

Introducing :

iTera - *Classic*



3 in 1 Core Technology

2021 4TH GENERATION ITERACARE DEVICE

An all new smarter, safer, more effective and portable focusing on partial and facial beauty therapy.

**“ Increases the Self Healing Ability in the Body
Activates Dormant Cells
Eliminates Unhealthy Cells
Clears Meridian and Lymph Points
Cleanses Blood impurities and Increases Blood Flow
Removes Body Humidity
Regulates Endocrine Glands and Protects our Organs**

When terahertz waves penetrate the body tissues they pass through healthy cells at a constant temperature. This generates millions of vibrations per second that can energize healthy cells down to the bone marrow which helps the body to self-heal and eliminates unhealthy cells.

Terahertz

QUANTUM



WHAT IS THE ITERACARE DEVICE

iTeraCare is based on energy healing and frequency. It addresses the illnesses caused due to disturbance of the energy flow in the body. It has passed all the safety tests – zero radiation, no electro-magnetic radiation and has human clinical reports. It is 100% safe, user friendly and can be used on anyone including new-born babies (for colic etc) and older people. You can use it to dry your hair rather than a hair dryer which has radiation and can cause headaches.



The iTeraCare Device uses:

Optical quartz crystal which draws out impurities from the body. It amplifies the energy and heat generated by the device and strengthens the auric field. It mimics 1/40 of sunlight.

Quantum energy Quantum is based on the principle of resonance and balance. Quantum generates scalar energy (healing energy) to balance the system and penetrates to the subatomic level.

Terahertz frequency resonates with our cells thus activating healthy and dormant cells while eliminating mutating and free radical cells in the process. This allows all our cells to function optimally. It is often referred to as “the wave of life”. Japanese scientists refer to it as the “light of life”

“If there’s a way to heat the bones, then all diseases can be treated”.

Ancient Greek Father of Medicine—Hippocrates

The iTeraCare wand uses targeted FREQUENCIES that can penetrate the body up to 20 - 30 centimetres, with zero discomfort, or irritation. Every cell in the body has a frequency measured in HERTZ. Terahertz resonates at the same FREQUENCY as healthy human, cells, generating millions of vibrations per second. The electromagnetic waves of the terahertz range are currently used in medicine as a biophysical factor for the correction of the several disorders.

FUNCTIONS:

Actively finds and protects the good cells which will increase energy;
Actively finds and eliminates mutated free radical cells & tumour cells;
Induces and strengthens the DNA organic molecules;
Clears meridian and lymph points, varicose veins, high blood pressure and is great for skin beautification (makes the skin firmer and clears skin disease);
Activates dormant stem cells in the bone marrow which will ultimately create new organs, tissues and muscles;
Improves microcirculation (in tiny nerves and arteries) – found in our brains and body that medicine and nutrition cannot reach;
Regulates the glands, endocrine glands and balances the hormones;
Removes body humidity – fats and water retention;
Protections against inflammation (the mother of all diseases).
Is excellent for relieving pain - muscles after workouts – shoulder pain, neck pain, lower back pain, headaches etc.
Regulates the body temperature and improves the immune system by activating weak and inactive cells.

The wave link can penetrate the skin between 20 to 30 cm which reaches the bone marrow, micro blood vessels & organs. It regulates the body temperature and improves the immune system by activating weak and inactive cells; It reduces blood flow resistance and viscosity.


Healthy blood cells promote smooth Chi and blood circulation, thus effectively transporting oxygen and nutrients to all body cells. Healthy Chi Meridian channels enhance organ recuperation, thus releasing body pain.

The device is similar to a hairdryer. There are 3 heat settings - you simply blow on acupuncture and meridian points and also on any areas of discomfort or pain. Before commencing treatment, you activate water by blowing it with the device. Drinking this charged water helps to eliminate toxins and flush the system during your treatment, it has shown to change the blood structure so the cells are separated, strong and moving fast. It can also kill free radicals. It has higher oxygen and restores the vitality of water. The water molecules are tiny and easily absorbed into your cells.

One hour of treatment is equivalent to 10 times of acupuncture; 20 times of cupping; and 40 times of massage.

The iTeraCare wand is the missing link for weight management because it easily helps one's body to increase metabolism, reduce water retention, increase circulation (blood and energy meridians), decrease inflammation, increase energy, increase nutrition bioavailability, relieve tension and stress, increase lymphatic drainage and detoxification and build immunity.

NOT RECOMMENDED FOR :



The iTeraCare device is suitable for:

Headaches; Sore throat/cough; colds; shoulder and neck pain, spleen, gastro; constipation; male prostate; athletes foot; varicose veins; rheumatic joint; high blood pressure; kidney issues, eyesight; hair problems; toothache; fever; bronchitis/asthma; lower back pain; female gynaecological; haemorrhoid; insomnia and sleep; hyperthyroidism; gout; obesity/weight loss; blood sugar imbalance; expelling cold and dampness; blockages in the meridian lines or lymphatic system; skin issues – acne, scars, blemishes, psoriasis; skin diseases; wrinkles (tightens the skin). Also can benefit Bell’s Palsy, Stroke, Alzheimer’s and Dementia.


The body becomes its own healer through the actions of the device.

iTeraCare  suitable for			
 Various Inflammations Rhinitis, Pharyngitis, Gastritis, Gynecological inflammation, Prostatitis, etc.	 Various Pains Shoulder and Neck pain, Joint pain, Low back pain, Headache, Dysmenorrhea, Sciatica, etc.	 Various Blockages Meridian blockage, Lymphatic blockage, Varicose veins, Weight Loss, Insomnia, Diabetes, High Blood Pressure etc.	 Various Skin Acne, Acne scars, Blemishes, Anti-aging, Beauty, Various Skin Diseases etc.

Terahertz Water can help improve the following:


- Acidity, gastric and intestinal problems; difficulty sleeping;
- Digestive problems—diarrhea, constipation, headaches;
- Liver problems—drowsiness, insomnia, itchy skin, abdominal pain;
- Kidney problems—low back pain, swelling of the legs, fatigue, increased blood pressure and lack of energy; Gynaecological diseases—unusual vaginal bleeding;
- Heart problems—rapid heartbeat, chest tightness, shoulder pain, dizziness, decreased blood pressure, anaemia and fever;
- Haemorrhoids—bleeding intestines; Hypertension—heart disease;
- Diabetes—frequent thirst, itchy skin, high blood sugar levels;
- Liver dysfunction—hepatitis, cirrhosis, itchy skin or rash, tinnitus.

According to animal studies, structured water can have an “increased rate of growth, reduced markers of oxidative stress, improved glycaemic and insulinemic responses in diabetics, improved blood lipid profile, improved semen and spermatozoa quality and increased tissue conductivity as measured using bioelectrical impedance analysis.”



Not only does terahertz frequency structure the water, but it can increase absorption of minerals by up to 50%.

The water energized by the terahertz device has a higher and more stable energy, can neutralize free radicals, increase oxygen content and restore the vitality of water. The iTeraCare device breaks the water particles into a linear form which allows it to be absorbed directly into our body cells.




Want to upregulate the "Guardian of your Genome" known as p53?

Studies show that the terahertz frequency, also called the "Wave of life", enhances cellular repair mechanisms by multiple proteins that play vital roles in DNA repair, including protein p53 (often called "a guardian of the genome"). That's right, terahertz frequency up-regulates the most studied gene in history!

The most commonly mutated gene in cancer, p53, works by shifting cell metabolism in a way that alters a cell's identity. More than 50% of all cancers contain a mutation in the p53 gene. The protein made from this gene is critical for ridding the body of cells that have developed genetic mistakes that could lead to cancer. Scientists have known about the p53 protein —

dubbed the "guardian of the genome" — since the 1970s. It is one of the most investigated molecules ever, the subject of more than 80,000 research papers.



- ❖ Cleanse blood impurities
- ❖ Cleanse inner body
- ❖ Promote blood circulation and removing clogs
- ❖ Soften blood vessels and lumps

STUDY—PEF AND SKIN REGENERATION

Skin is the largest organ in the body and has important functional and psychological significance. Aging, trauma, and chronic metabolic diseases, including diabetes, often lead to alterations in skin colour, texture, and barrier function. Loss of skin function leads to atrophy, infection, chronic wounds, laxity, and rhytids. Epidermal atrophy – characterized by a thinning of the epidermis and an increase in fragility – is observed in 32% of individuals over 60 years of age. Thus, it is not surprising that 2,156,075 skin rejuvenation procedures were performed in 2013 in the United States alone.

Because these skin diseases are often associated with poor re-epithelialization, poor blood supply, reduced collagenases and a loss of collagen functional properties, current rejuvenation therapies are focused on the removal of non-functional tissue and the induction of de novo growth of healthy dermal cells, blood vessels, and extracellular matrix.

Therapies that aim to induce collagen synthesis and to restore skin function and appearance include non-ablative laser technologies, radiofrequency, ultrasound, electro-optical synergy, chemical peels, microdermabrasion, injectable fillers, neurotoxins, skin needling, mesotherapy, platelet-rich plasma, and cell therapies.



And now Quantum Terahertz Frequency Physiotherapy (pulsed electric fields). Researchers report that certain PEF regimens applied non-invasively lead to epithelial growth, an increase in vascular supply to the skin, the formation of new microvasculature, and increased collagen deposition without scarring.

iTeraCare Principles of Therapy 5 in 1

The image shows two handheld devices: a white iTera-Pro on the left and a red iTera-Classic on the right. Between them is a list of five benefits in colored boxes. To the right of the iTera-Classic is a text box comparing the therapy to other treatments.

- Activate and Repair Cells
- Rapidly Clear Meridians
- Rapidly Detoxify Lymphatic System
- Rapidly Remove Water Retention
- Rapidly Improve Blood Microcirculation

An hour of therapy is equal to:
10 times acupuncture
20 times of cupping
30 times of Gua Sha
40 times of massage

iTera-Pro iTera-Classic

HOW TO OPTIMISE RESULTS

1. Drink 2 glasses of charged warm water before and after the therapy. For those who can't drink the recommended amount, you may drink in between or drink staggered within an hour or 2 during and after the session.
2. For the first 3 therapies follow the 10 step warm up guide to slowly introduce the THZ Frequency to your body. This will clear blockages and open meridian points. This is a total of 20 mins per 10 step guide session.
3. Next, focus on uncomfortable areas, or areas you want to treat and activate cellular generation and elimination of bad cells. 5-15 mins. per area.
4. Never blow on areas where there's haemorrhage, open wounds, cuts or burns; areas where there's cancer cells nearby blood vessel causing bleeding. This may cause more pain and unwanted reactions.
5. Continue taking your supplements and medicine while using the iTeraCare device. Once you see significant progress you can cut down taking your meds provided you consult your Doctor.
6. Charging your tea, natural fruit juice etc. is perfect to increase detoxification effect and nutrients in the body. Never charge carbonated drinks or water in plastic container. You can also charge your food.
7. We recommend a total of 15-30 minutes per day for healthy individuals. Not more than 2 hours per day for those having health challenges. Take a break for up to 3 days every month or 1 day per week after continuous use. The break is important so the body's cells won't get too immune with the frequency.
8. Once you undergo a healing crisis, we recommend changing your lifestyle to a healthier habit. Take enough rest, drink up to 4 litres of water per day, have proper nutrition, exercise, sunlight and deep breathing. Cut down the use of the device up to 50% of usual time, eg. from 30 mins. down to 15 mins. until you feel well again.
9. Once you achieve your health goal, minimize regular prolonged use of the device. Using it for 5-10 mins per day is perfect. Take a break once in a while. But don't stop drinking lots of charged water.
10. Removing your clothes during the therapy is highly recommended to maximize the iTeraCare heat effect in the body.
11. Never eat raw meat and cold foods within 4 hours after the therapy.
12. Don't take a bath within 4 hours after the therapy. No electric fan and air-condition directly pointed on you within 4 hours after the therapy also.
13. Enjoy the journey of healing and optimal health. Never worry and always be joyful for within you is the power to heal yourself.

TIPS:

1. Try to blow it yourself—because you can feel and control the penetration.
2. Try to blow directly on the skin—its also okay to blow through clothing.
3. You must be patient for your body to repair itself, Let the pain you have accumulated for years, slowly be blown away.
4. When blood cells are healthy, blood circulation improves. It can now effectively bring oxygen and nutrition to other cells and flush away toxins and dead cells. Then a lot of pain can be eliminated.
5. When the meridian cells are healthy, the nervous system will be unblocked and the organs can be nourished to function effectively.
6. Every time you blow, you should blow on the 3 acupuncture points, the 5 lymphatic nodes, and the 2 main meridians, as well as the any sore areas.
7. It is best not to get into cold air-conditioning for a few hours after therapy, or to eat cold food or drinks. Also do not take a shower.
8. Be sure to understand the contraindications and precautions before use.
- 9/ Disclaimer: The above suggestions are only personal opinions—please make your own judgment.



What is the difference between the iTera Classic & iTera Pro?

The Classic is for personal use. It needs to cool down between uses and is not suitable for use all day as would occur in a professional setting. It penetrates 7-11 inches, 50 hz, 800 watts and 3 modes (low, high and high with heat), runs for 45 minutes without a break.

The PRO wand is slightly bigger and much stronger. It works faster and more softly. You don't need to let it cool down between uses. It penetrates 7-11 inches, 50 hz, 1200 watts, 4 modes (low heat, high heat, low cool, and high cool), quieter model, Cool is great for facial applications, runs for 24 hours without a break and treatment times are lower. Perfect for professionals.


IteraCare is also next level when it comes to the spiritual / meditation possibilities. In essence it aligns you with the frequency of the earth - 7.6THz - and can put you into an Alpha state. Both machines do this.

Steps of using iTeraCare™

Warm up the body. Blow step 1 to 10, each part for about 2-3mins



Drink 2 glasses
of warm water



1. Palms




2. Soles



3. Behind
the Knees



4. Groin



5. Armpit



6. Lymph




7. Upper Back



8. Head



9. Butt



10. Blow Downwards




Drink 2 more glasses
of warm water

iTeraCare device can also be blown on Head, Face, Body, Limbs and Pain areas.

*High Blood Pressure user don't blow no. 8



iTeraCare Heat



BENEFITS OF HEATING THE BODY


- Boost Immune System**
- Increases Metabolic Activity in Cells**
- Increases Cellular Oxygenation**
- Blood Flow Increases**
- Metabolizes Waste Products in the Body**
- Softens Connective Tissues**
- Relaxes Muscle Fibers**
- Reduces Joint Stiffness**
- Powerful Pain Killer**

CELLULITE DETOX PROTOCOL


Cellulite is often thought of as unsightly, but it can also be indicative of other health issues. Often times fat deposits are created when the kidneys, liver and digestive system fail to adequately eliminate waste products. Adding this protocol daily might offer some benefits. Regarding the duration, start slow and increase as your body tolerates. The heat setting is ideal, but you may have to work up to using high heat. Please consult your medical provider if you have any questions. This is not medical advice, but several people have had amazing results.

Protocol Do the basic protocol; Drink plenty of charged water throughout the day; Charge your throat to help regulate the hormones as estrogen has a large role in cellulite formation.; Charge your stomach in a clockwise pattern

to help open the digestive pathways for elimination. Make sure you are eliminating daily to release toxins; Lifestyle changes; clean diet, eliminate hormone disrupting foods, get rid of plastics and other toxins in your environment, stay hydrated, exercise daily.



Vertigo and Dizziness Protocol



This protocol seems to diminish all of the symptoms and they gradually dissipate.

Protocol: Do the basic protocol; Charge the points above; Repeat daily. Start with 30 seconds on each spot to see how your body tolerates the Terahertz frequency and build up to 1-2 min on each spot. You can also incorporate the spine.

Contraindications: do not blow on the top of the head (GV 20) if someone has high blood pressure. Do not blow on metal or implants. It is always best to make sure that these symptoms are not part of a larger problem. Please consult a physician if needed or symptoms persist.

This is not medical advice.

Terahertz waves vibrate at the same frequency as normal body cells. When terahertz waves penetrate the body tissues, they pass through healthy cells at constant temperature, however, they cause unhealthy cells to vibrate faster which raise their cellular temperature unusually higher, thus destroying these unhealthy cells in the process. Terahertz waves together with quantum technology, can penetrate deeper and faster.

-TCM Wellness

TESTIMONIALS



SUSAN RENNIE Lately my stomach has been bubbling and burping, and it always felt I had to hold it. After 1 session with the wand, my stomach felt “calm” and that it had again become part of my body. No more burping either. I had knocked my finger whilst moving a trailer and the knuckle was swollen and very painful. The following morning my finger was pain free and the swelling had gone away. The iTeraCare water seems to have a different consistency and my cells feel like they have been hydrated. All this from one session. Cant wait to get mine.


WENDY SCOTT. I am a Wellness Therapist and help to relieve pain and emotion blockages from clients.. I have been using the wand for the last couple of weeks and it has changed my life for the better. No more brain fog, I can think clearly now and my memory is starting to come back. I have also noticed that my spider veins on my legs are disappearing, and skin feels amazing. I recommend this wonder machine to all my clients.



SANDRA BARNSELEY I have been painting the outside of my house and I moved a ladder not realizing there was a heavy battery and drill sitting on top of it. It fell from a height above my head and landed on my big toe. The pain was excruciating, and I thought I had broken my toe. I went straight inside and did the iTeraCare on it for a few minutes and the pain went away and there has been no pain at all from that moment.

I also had a mole that was painful if I touched it. I have tried all sorts of potions with no change. I used iTeracare on it for 30 seconds and the pain went. Gone, nada, zip nothing. That was what really made me realize what I had in my hands.

SARA HOME Izzy kindly did an iTeraCare session on me for my shingles on day 4. The affected area was very red and angry looking. It was sore and I had added nerve pain in areas of my skin not directly having blisters. The day after the session was tough with detoxing. The following day the redness had begun to subside. I'm writing this on day 6th. Far less nerve pain and redness, and much more energy. I'm absolutely convinced iTeraCare is helping my body move through this much much quicker. Much gratitude to Izzy and Marita





DAVID CLAYTON A rapidly expanding growth in the corner of my eye, the size of a pea was treated twice with the iTeraCare wand over 48 hours. I then had a portion removed by a doctor to get tested. A week later, I had a third session ... the next day I was overjoyed to discover the growth had disappeared completely. You can see why I highly recommend this device.

PAULA MACK When I first used the device I had hurt my big toe which was very sore and black underneath. I could barely touch it or the nail. Also, my whole body was very sore as I was on day 4 of a 30 day exercise challenge, so all my muscles were achy. After a 20 minute session with the device I was so amazed that my toe and nail were no longer sore. I could actually touch them. Also the pain in my body had also gone. Afterwards I went to the supermarket and was walking around feeling like I was in my 20's again and that feeling of being so in love with life. It was a very profound experience! I'm so excited about this device.



MARLENE GALLER I have been using the iTeraCare wand for the past 4 weeks. I have had an arthritic knee for 10 years and the pain disappeared instantly and then within a week I was able to kneel on the floor with my grand daughter with no pain at all. My hip and lower back become more flexible and stable each day. I can now go for walks without pain which I haven't been able to do for the last 4 years. I have more energy everyday, I sleep deeper and need less food and my body shape continually keeps changing for the better. Most of all, I simply love being able to feel the deep relaxation it brings to my whole being. It feels amazing to be in control of my own body's health and vitality

MARITA MASON I have been using the wand for 5 weeks and am feeling sensational. I haven't been able to play tennis for a year because of a chip in my knee and I had been very sluggish. I'm prone to constipation but that changed over night and I now don't get headaches from drinking red wine and coffee or eating chocolate. The chip in the knee kept me awake at night but that is now better. I need 2 hours less sleep and wake up super refreshed. My skin is improving, cellulite is dissolving, body is toning up and my nails are getting stronger. My "chicken neck" and "granny wings" are also tightening up and I have so much more life force, enthusiasm and motivation. Who would have thought you could feel so great in your 70's? I'm so very grateful.





LISA MURRAY Rotator cuff pain gone. Throat clearing greatly diminished. Skin noticeably improving. Big energetic transformations. In just over two weeks of intermittent use, I've noticed all of these things personally. I've also done sessions on 30+ people and they have reported instant pain relief, a huge feeling of space and incredible meditation experiences. I can't wait for my wand to arrive!

John Scevola The iTeraCare device has definitely helped me personally to improve my weight management. It has helped lower my blood pressure, eliminate acid reflux, eliminate lazy leg syndrome, improve my sleeping through the night, eliminate excruciating lower back discomfort, and eliminate a tooth abscess.



2 months before the therapy



OSTEOPOROSIS

93 year old lady with severe osteoporosis.

Nerves impingement causing great pain. Unable to even move an inch. She has to remain in bed. This was 2 months ago. Now today, she can even do the track cycle after doing our iteracare blower therapy twice a day for 2 months till today. She can even lift her leg to do karaoke with great ease. Marvelous recovery. Thanks to our wonderful iteracare blower.



Friendly Reminder

Reaction :
Using iTeraCare Therapy device may have some reactions:

1. Due to different individual physical body structure, some may experience painful on certain areas.
2. Previous unhealed or untreated pains or symptoms will appear during treatment.
3. During therapy, such reaction may appear from time to time.


All of the above reactions are normal due to its therapy healing process.

*** After Therapy: For fast recovery and healing, you are strongly advice NOT to consume cold water and/ or cold food, or eat raw, No Fan, No Air Conditioner. You can only shower after 4 hours of the therapy, otherwise it will be counterproductive!



SLE (Lupus) After 2 sessions of therapy






Ms. Aida Patana from Philippines. She got severe Covid-19, admitted 25days in the hospital. She was breathing thru the help of oxygen tank. According to the doctors, her recovery might take 3 months due to severe complications. 😞 Then she discovered iTeraCare Device. After her 1st 40-minute therapy, she was able to breathe comfortably and used it for 1-week straight. now she's finally well, no more oxygen tank and took her less than 2 weeks to fully recover from Covid-19. Thank you Iteracare Device. 🙏👉👈



Breast Cancer 30% lump reduced; 90% pain relieved within 22 days of therapy.



**10 years Vasculitis—
solved after
iTeraCare sessions**



Iteracare


After losing weight there was a lot of sagging skin. After treatment on his left leg, the skin became taut.

**Before
iteracare** 🙏




**After 5 minutes
ITeraCare** 🙏

Psoriasis



After 7 days



After 9 days



After 22 days




PARALYSED PATIENT started walking after 15 treatments. Improvement in circulation in 5 days.



Bell's Palsy after a week of therapy



SEVERE BACK PAIN. After 3 days of treatment the pain had totally gone.



The Advantage of iTeraCare^W

1.	Manufacturer is continuously focusing and upgrading the 4 th Generation terahertz technology
2	Based on 3 in 1 Core Technology 1. Terahertz technology 2. Quantum technology 3. Optical Quartz Technology
3	Passed strict inspection by the National Authority for its safe and efficacy and effectiveness
4	Numerous proven test reports from the Chinese Academy of Sciences for the physical therapy effectiveness
5	Quality Assurance of One year warranty
6	National Renowned, Certified and Patented
7	iNetwork is the Sole Distributor for Ten countries worldwide

GENERAL DISCLAIMER

The iTeraCare device is a certified household electrical device. It has been certified as “safe to use”. It emits no harmful radiation and it has labels allowing for its distribution internationally.

The iTeraCare blower is not a medical device and the manufacturer provides no medical advice. Users understand that the device is not intended to be used as a medical device or for diagnostic purposes on which reliance should be placed. Its use is not intended to be relied upon in lieu of medical treatment or advice by a trained medical care provider. Users understand that they use the iTeraCare blower at their own risk.

Always consult your doctor or other healthcare professional with any questions regarding any medical condition, or for specific guidance regarding the optimal use of the iTeraCare blower, nutrition or physical activity.

I believe this device has been brought to the world at a time when it is most needed. There is so much confusion, division and ill health at the moment. For 5 weeks I've seen so many peoples' lives transformed—whether its emotionally, spiritually or physically through the simple process of “blowing”. The device is affordable for most and can be used on every member in the household—even babies and the elderly, with no special skills required. As people heal, I believe we will come together as a loving collective without separation and fear. I see this as the beginning of a new era on earth.

Marita Mason
HEALING of the HEART




Sd truth



The future of medicine is **FREQUENCY**

Do you know that our soul is composed of harmony?
You can look at disease as a form of disharmony.
And there's no organ system in the body that's not
affected by sound, music and vibration.



Future medicine will
be the medicine of
frequencies.

Albert Einstein

If you want to supercharge
your health with this “magic
wand” PLEASE CONTACT: